



PRODUCT SPECIFICATIONS

Product Name:	PG II protein	
Description:	Recombinant Human Pepsinogen II protein with a C-terminal His tag	
Catalog #:	A020602R	Lot #: AY22M2001
Source:	<i>HEK293 Cells</i>	
Format:	Purified, Liquid	
Concentration:	1.0mg/mL (+/- 10 %)	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	10%Glycerol, pH7.4.	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular application.	
Storage:	Store at 2-8°C for two weeks and -20°C to -80°C for long time. Avoid multiple freeze/thaw cycles.	
production date:	01/2022	
Expiry Date:	01/2025	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com