



PRODUCT SPECIFICATIONS

| | | |
|-------------------------|---|-------------------------|
| Product Name: | MAb to PIVKA-II | |
| Description : | Monoclonal Antibody to Protein Induced by Vitamin K Absence or Antagonist-II (PIVKA-II) | |
| Catalog #: | M020902M | Lot #: MY22M3003 |
| Isotype: | IgG1 | |
| Immunogen: | Purified Protein Induced by Vitamin K Absence or Antagonist-II (PIVKA-II) | |
| Source: | Mouse ascites | |
| Format: | Purified, Liquid | |
| Concentration: | 4.51 mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4) | |
| Purity: | ≥90% (SDS-PAGE) | |
| Product Buffer: | PBS, pH 7.4 | |
| Preservative: | 0.05% Proclin 300 | |
| Applications: | Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications. | |
| Storage: | Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles. | |
| Production Date: | 03/2022 | |
| Expiry Date: | 03/2025 | |

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu
Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com