



PRODUCT SPECIFICATIONS

Product Name:	MAb to PIVKA-II	
Description :	Monoclonal Antibody to Protein Induced by Vitamin K Absence or Antagonist-II (PIVKA-II)	
Catalog #:	M021104M	Lot #: MY22M1102
Isotype:	IgG1	
Immunogen:	Purified Protein Induced by Vitamin K Absence or Antagonist-II (PIVKA-II)	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration:	2.26 mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	PBS, pH 7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
Production Date:	02/2022	
Expiry Date:	02/2025	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu
Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com