



PRODUCT SPECIFICATIONS

Product Name:	MAb to ProGRP	
Description:	Monoclonal Antibody to Pro-gastrin-releasing peptide	
Catalog #:	M021902M	Lot #: [Lot specific]
Isotype:	IgG1	
Immunogen:	Purified ProGRP protein	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration :	[Lot specific] mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity :	≥90% (SDS-PAGE)	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
	Capture	Detection
	M021902M	M021901M/M021904M/M021905M/M021906M
Storage:	Store at -20°C . Avoid multiple freeze/thaw cycles.	
Expiry Date:	[Lot specific]	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com