



PRODUCT SPECIFICATIONS

Product Name:	MAb to SCCA	
Description:	Monoclonal Antibody to Squamous Cell Carcinoma Antigen	
Catalog #:	M021210M	Lot #: MY22M1304
Isotype:	IgG1	
Immunogen:	Purified SCCA	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration :	1.7mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity :	≥90% (SDS-PAGE)	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	

Capture

M021205M

Detection

M021206M/M021207M/M021208M/

M021209M/M021210M/M021211M

Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.
-----------------	---------------------------------------------------------------------------------------------

Production Date: 04/2022

Expiry Date: 04/2025

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com