




PRODUCT SPECIFICATIONS

Product Name:	VEGF-165 Protein	
Description:	Vascular endothelial growth factor Antigen with C-terminal HIS tag	
Catalog #:	A051701R	Lot#: AY22M1502
Source:	<i>Escherichia coli</i>	
Format:	Purified, Liquid	
Concentration:	1.0 mg/mL(+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	PBS, pH7.4	
Preservative:	0.05% ProClin300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular application.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
Production Date	02/2022	
Expiry Date:	02/2025	

SDS-PAGE:



Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com