



PRODUCT SPECIFICATIONS

Product Name:	ALD-HRP Conjugate	
Description:	Aldosterone HRP Conjugate; Aldosterone Antigen	
Catalog #:	A051002R	Lot #: AY21M0807
Source:	Synthesis	
Format:	HRP, Liquid	
Concentration:	7.0 mg/mL (+/- 10 %)	
Product Buffer:	10mM PBS, pH7.4, 50% glycerol,	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
Production Date:	07/2021	
Expiry Date:	07/2024	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com