



PRODUCT SPECIFICATIONS

Product Name:	MAb to Cortisol	
Description :	Monoclonal Antibody to Hydrocortisone	
Catalog #:	M051901M	Lot #: [Lot specific]
Isotype:	IgG2a	
Immunogen:	Purified Hydrocortisone	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration:	[Lot specific] mg/mL (+/- 10 %) (OD280nm, $E^{0.1\%} = 1.4$)	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	PBS , pH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
Production Date:	[Lot specific]	
Expiry Date:	[Lot specific]	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com