



PRODUCT SPECIFICATIONS

| | | |
|-------------------------|---|-------------------------|
| Product Name: | VPA-BCP Conjugate | |
| Description : | Valproic Acid BCP Conjugate | |
| Catalog #: | A151701R | Lot #: AY21M1704 |
| Source: | Synthesis | |
| Format: | Liquid | |
| Concentration: | 5 mg/mL (+/- 10 %) | |
| Purity: | None | |
| Product Buffer: | PBS, pH7.4 | |
| Preservative: | 0.05% Proclin 300 | |
| Applications: | Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications. | |
| Storage: | Store at 2-8°C. | |
| production date: | 04/2021 | |
| Expiry Date: | 04/2024 | |

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com