



PRODUCT SPECIFICATIONS

Product Name:	MAb to FCoV	
Description:	Monoclonal Antibody to Feline Coronavirus	
Catalog #:	M101203M	Lot #: MY21M0412
Isotype:	IgG	
Immunogen:	Purified FCoV protein	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration:	12.3mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity:	≥95% (SDS-PAGE)	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C	
Expiry Date:	12/2024	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY