



PRODUCT SPECIFICATIONS

Product Name:	CDV N Protein	
Description:	Canine Distemper Virus N Protein with C-terminal HIS tag	
Catalog #:	A100302R	Lot #: AY22M2508
Source:	<i>Escherichia coli</i>	
Format:	Purified, Liquid	
Concentration:	1.5mg/mL (+/- 10 %)	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	Tris , PH8.5	
Preservative:	0.05%Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time . Avoid multiple freeze/thaw cycles.	
Production date:	08/2022	
Expiry Date:	08/2025	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com