



## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	MAb to CDV	
<b>Description:</b>	Monoclonal Antibody to Canine Distemper Virus	
<b>Catalog #:</b>	M100302M	<b>Lot #:</b> MY21M0709
<b>Isotype:</b>	IgG1	
<b>Immunogen:</b>	Purified CDV protein	
<b>Source:</b>	Mouse ascites	
<b>Format:</b>	Purified, Liquid	
<b>Concentration :</b>	8.10 mg/mL (+/- 10 %) (OD280nm, E <sup>0.1%</sup> = 1.4)	
<b>Purity :</b>	≥90% (SDS-PAGE)	
<b>Product Buffer:</b>	PBS, PH7.4	
<b>Preservative:</b>	0.05% Proclin 300	
<b>Applications:</b>	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
<b>Storage:</b>	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
<b>Production Date:</b>	09/2021	
<b>Expiry Date:</b>	09/2024	

Web: [www.firsunbioscience.com](http://www.firsunbioscience.com)

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003      Email: [order@firsunbioscience.com](mailto:order@firsunbioscience.com)



**Note:** Centrifuge before opening to ensure complete recovery of vial contents

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

Web: [www.firsunbioscience.com](http://www.firsunbioscience.com)

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003      Email: [order@firsunbioscience.com](mailto:order@firsunbioscience.com)