



## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	MAb to CPV	
<b>Description:</b>	Monoclonal Antibody to Canine Parvovirus	
<b>Catalog #:</b>	M100207M	<b>Lot #:</b> MY22M0702
<b>Isotype:</b>	IgG	
<b>Immunogen:</b>	Purified CPV protein	
<b>Source:</b>	Mouse ascites	
<b>Format:</b>	Purified, Liquid	
<b>Concentration :</b>	2.68 mg/mL (+/- 10 %) (OD280nm, E <sup>0.1%</sup> = 1.4)	
<b>Purity :</b>	≥90% (SDS-PAGE)	
<b>Product Buffer:</b>	PBS, PH7.4	
<b>Preservative:</b>	0.05% Proclin 300	
<b>Applications:</b>	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
<b>Storage:</b>	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
<b>production date:</b>	02/2022	
<b>Expiry Date:</b>	02/2025	

**Note:** Centrifuge before opening to ensure complete recovery of vial contents

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

Web: [www.firsunbioscience.com](http://www.firsunbioscience.com)

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003      Email: [order@firsunbioscience.com](mailto:order@firsunbioscience.com)