



PRODUCT SPECIFICATIONS

Product Name:	RV Protein	
Description:	Rabies Virus Protein	
Catalog #:	A010502N	Lot #: AY23M0103
Source:	<i>BHK Cell</i>	
Format:	Purified, Liquid	
Concentration:	2.7 mg/mL (+/- 10 %)	
Purity:	None	
Product Buffer:	PBS, PH7.4	
Preservative:	1‰ Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Inactivation:	Whole infective virus is inactivated with formalin. Confirm by attempting growth under original culture conditions.	
Storage:	Store at -80°C. Avoid multiple freeze/thaw cycles.	
Production Date:	03/2023	
Expiry Date:	03/2026	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com