



## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	Recombinant Canine CRP protein	
<b>Description:</b>	Recombinant Canis CRP protein with C-terminus HIS tag	
<b>Catalog #:</b>	A100501R	<b>Lot #:</b> AY22M2001
<b>Source:</b>	HEK 293 Cells	
<b>Format:</b>	Purified, Liquid	
<b>Concentration:</b>	1.0 mg/mL (+/- 10 %)	
<b>Purity:</b>	≥90% (SDS-PAGE)	
<b>Product Buffer:</b>	PBS, pH 7.4, 3%trehalose ,0.01%Tween-80	
<b>Applications:</b>	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
<b>Storage:</b>	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
<b>production date:</b>	01/2022	
<b>Expiry Date:</b>	01/2025	

**Note:** Centrifuge before opening to ensure complete recovery of vial contents

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

Web: [www.firsunbioscience.com](http://www.firsunbioscience.com)

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003      Email: [order@firsunbioscience.com](mailto:order@firsunbioscience.com)