



PRODUCT SPECIFICATIONS

Product Name:	MAb to Canine CRP	
Description:	Monoclonal Antibody to Canine C-Reactive Protein	
Catalog #:	M100504M	Lot #: MY21M1712
Isotype:	IgG1	
Immunogen:	Purified Canine C-Reactive Protein	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration :	2.36 mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity :	≥90% (SDS-PAGE)	
Product Buffer:	PBS, pH 7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
production date:	12/2021	
Expiry Date:	12/2024	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com