



PRODUCT SPECIFICATIONS

Product Name:	MAb to ASFV (P72)	
Description :	Monoclonal Antibody to p72 Protein of African Swine Fever Virus	
Catalog #:	M010804M	Lot #: MY21M0909
Isotype:	IgG1	
Immunogen:	Purified recombinant VP72 protein	
Source:	Mouse ascites	
Format:	Purified, Liquid	
Concentration :	5.05 mg/mL (+/- 10 %) (OD280nm, E ^{0.1%} = 1.4)	
Purity :	≥90% (SDS-PAGE)	
Product Buffer:	PBS, pH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular application.	

Capture	Detection
M010801M	M010804M

Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.
Production Date:	09/2021
Expiry Date:	09/2024

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com