




## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	CCV N Protein	
<b>Description:</b>	Canine Coronavirus N Protein with C-terminal HIS tag	
<b>Catalog #:</b>	A100401R	<b>Lot #:</b> AY23M2605
<b>Source:</b>	<i>Escherichia coli</i>	
<b>Format:</b>	Purified, Liquid	
<b>Concentration:</b>	1.5 mg/mL (+/- 10 %)	
<b>Purity:</b>	≥90% (SDS-PAGE)	
<b>Product Buffer:</b>	Tris PH7.5	
<b>Preservative:</b>	1‰ Proclin 300	
<b>Applications:</b>	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
<b>Storage:</b>	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.	
<b>production date:</b>	05/2023	
<b>Expiry Date:</b>	05/2026	
<b>SDS-PAGE:</b>		

**Note:** Centrifuge before opening to ensure complete recovery of vial contents

### FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: [www.firsunbioscience.com](http://www.firsunbioscience.com)

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003      Email: [order@firsunbioscience.com](mailto:order@firsunbioscience.com)