



PRODUCT SPECIFICATIONS

Product Name:	PAb to Rabbit anti-Cat IgG	
Description:	Rabbit Polyclonal Antibody to Cat Immunoglobulin G (H+L)	
Catalog #:	T111601R	Lot #: PY22M0207
Immunogen:	Purified Cat IgG1	
Source:	Rabbit Serum	
Format:	Purified, Liquid	
Concentration:	3.70 mg/mL	
Purity:	≥90% (SDS-PAGE)	
Product Buffer:	PBS, PH7.4	
Preservative:	0.05% Proclin 300	
Applications:	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
Storage:	Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles. Avoid Light.	
Production Date:	07/2022	
Expiry Date:	07/2025	

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com