



PRODUCT SPECIFICATIONS

| | | |
|-------------------------|---|-------------------------|
| Product Name: | SDMA-Biotin Conjugate | |
| Description: | Symmetric dimethylarginine Conjugate with Biotin | |
| Catalog #: | A100902R | Lot #: AY23M1410 |
| Source: | Synthesis | |
| Format: | Biotin, Liquid | |
| Concentration: | 5.1 mg/mL (+/- 10 %) | |
| Purity: | N/A | |
| Product Buffer: | PBS, pH7.4 | |
| Preservative: | 0.05% Proclin 300 | |
| Applications: | Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications. | |
| Storage: | Store at 2-8°C in dark. | |
| Production Date: | 10/2023 | |
| Expiry Date: | 10/2026 | |

Note: Centrifuge before opening to ensure complete recovery of vial contents

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Web: www.firsunbioscience.com

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003 Email: order@firsunbioscience.com