



## PRODUCT SPECIFICATIONS

<b>Product Name:</b>	SDMA-OVA Conjugate	
<b>Description:</b>	Symmetric dimethylarginine Conjugate with OVA	
<b>Catalog #:</b>	A100903R	<b>Lot #:</b> AY23M2407
<b>Source:</b>	Synthesis	
<b>Format:</b>	OVA, Liquid	
<b>Concentration:</b>	1.6 mg/mL (+/- 10 %)	
<b>Purity:</b>	N/A	
<b>Product Buffer:</b>	PBS, pH7.4	
<b>Preservative:</b>	0.05% Proclin 300	
<b>Applications:</b>	Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications.	
<b>Storage:</b>	Store at 2-8°C in dark.	
<b>Production Date:</b>	07/2023	
<b>Expiry Date:</b>	07/2026	

**Note:** Centrifuge before opening to ensure complete recovery of vial contents

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

Web: [www.firsunbioscience.com](http://www.firsunbioscience.com)

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003      Email: [order@firsunbioscience.com](mailto:order@firsunbioscience.com)