



## PRODUCT SPECIFICATIONS

|                         |   |                         |
|-------------------------|---|-------------------------|
| <b>Product Name:</b>    | MAb to RBC  |                         |
| <b>Description:</b>     | Monoclonal Antibody to Red Blood Cell   |                         |
| <b>Catalog #:</b>       | T112601M  | <b>Lot #:</b> MY21M3012 |
| <b>Isotype:</b>         | IgG1  |                         |
| <b>Immunogen:</b>       | Purified Human RBC  |                         |
| <b>Source:</b>          | Mouse ascites   |                         |
| <b>Format:</b>          | Purified, Liquid  |                         |
| <b>Concentration:</b>   | 10mg/mL (+/- 10 %) (OD280nm, E <sup>0.1%</sup> = 1.4)   |                         |
| <b>Purity:</b>          | ≥90% (SDS-PAGE)   |                         |
| <b>Product Buffer:</b>  | PBS, PH7.4  |                         |
| <b>Preservative:</b>    | 0.05% ProClin300  |                         |
| <b>Applications:</b>    | Suitable for use in immunoassays. Each laboratory should determine an optimum working titer for use in its particular applications. |                         |
| <b>Storage:</b>         | Store at 2-8°C for two weeks and -20°C for long time. Avoid multiple freeze/thaw cycles.  |                         |
| <b>Production Date:</b> | 12/2021   |                         |
| <b>Expiry Date:</b>     | 12/2024   |                         |

**Note:** Centrifuge before opening to ensure complete recovery of vial contents

### **FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

Web: [www.firsunbioscience.com](http://www.firsunbioscience.com)

Add: Building C5, Hongfeng Science and Technology Park, Qixia District, Nanjing, Jiangsu Province, China, 210000

Tel:86-025-83108003      Email: [order@firsunbioscience.com](mailto:order@firsunbioscience.com)